

AL-02



AL-02 contains AP-9925, ZFB-01.

AL SERIAL SOLVENT FILTRATION APPARATUS

Widely used in gravimeter, microanalysis, colloid dispersion and bacteria-free experiment, the AL-01 Solvent Filtration Apparatus filtrates and degas the mobile phase of HPLC analysis. Playing a critical role in lifetime dilation of the apparatus and HPLC column, the AL-01 Solvent Filtration Apparatus can enhance measurement accuracy.

Filtration Bottles

These elegant Filtration Bottles comes with many advantages as a result of their high-quality and hard-glass materials. All filtration bottles have high-pressured endurance, fabulous interchange ability, and well-proportional wall. They come in international standard sizes that match most of the famous foreign brands of products.

The Vacuum/Pressure Pump

The Vacuum Pressure Pump has high-working efficiency while producing very little noise. Offered in two types of vacuum: Steady Pressure Pump and Adjustable Pressure Pump.

Solvent Pickup Adaptor

The newly designed Solvent Pickup Adaptor filtrates and/or collects solvent directly from solvent bottle through film of filter. This new design allows the process of filtrating easier, dramatically increases effectiveness, and reduces air pollution.