

## AP-01D



Type :	AP-01D
Max vacuum (Kpa):	0-80 Adjustable
Max Pressure(psi):	_____
Free-air Capability (l/min):	8
Power (VA):	50
Power supply:	AC220V/50Hz
Dimension (mm) :	246Lx236Wx136H
Shpg (Kg) :	3.1

### Advantage

Main parts are made of Teflon, PTFE resin.

Can endure large volume flow and high pressure and meet common filter systems' requirements.

Have positive and negative connectors separately and can be used in many areas.