

AP-01P



Type :	AP-01P
Max vacuum (Kpa): :	≥80KPa
Max Pressure(psi):	_____
Free-air Capability (l/min):	8
Power (VA):	50
Power supply:	AC220V/50Hz
Dimension (mm) :	246Lx190Wx135H
Shpg (Kg) :	3.2 Advantage

Main parts are made of Teflon, PTFE resin.

Can endure large volume flow and high pressure and meet common filter systems' requirements.
Have positive and negative connectors separately and can be used in many areas.