

## AP-02B



Type:	AP-02B
Max vacuum (Kpa):	0-80 Adjustable
Max Pressure(psi):	18
Free-air Capability (l/min):	10
Power (VA):	45
Power supply:	AC220V/50Hz
Dimension (mm) :	230Lx85Wx200H
Shpg (Kg) :	2.5

The newest product from Auto science, the AP-02B is an adjustable pressure pump with a unique method of adjusting pressure. Using the pressure sender to detect the pressure in the system and then re-adjust the pump's working period, the AP-02B maintains a stable pressure while the current changes. With adjustable stable pressure and long-time working, the AP-02B is especially suitable for working condition as it works intermittently majority of the time.

- Its flat shape body is perfect for placement on lab benches.
- Accurate, reliable dual 'LED' pressure displays.
- Convenient knob for adjustment of vacuum pressure.
- Positive and negative joints allow the system to be widely applicator.
- Double the current through its dual pump-head design (AP-02B).