

AP-9901S



Type:	AP-9901S
Max vacuum (Kpa):	≥ 70 KPa
Max Pressure(Kpa):	> 150
Free-air Capability (l/min):	12
Power (VA):	80
Power supply:	AC220V/50Hz
Dimension (mm) :	210Lx140Wx120H
Shpg (Kg) :	3.5Advantage

Main parts are made of Teflon, PTFE resin.

Can endure large volume flow and high pressure and meet common filter systems' requirements.

Have positive and negative connectors separately and can be used in many areas.