

## AP-9950



Type :	AP-9950
Max vacuum (Kpa):	0-80 Adjustable
Max Pressure(psi):	20
Free-air Capability (l/min):	50
Power (VA):	180
Power supply:	AC220V/50Hz
Dimension (mm) :	269Lx210Wx155H
Shpg (Kg) :	8.2

### Advantage

Main parts are made of Teflon, PTFE resin.

Can endure large volume flow and high pressure and meet common filter systems' requirements.

Have positive and negative connectors separately and can be used in many areas.