

## API-1 Autosampler



### The specification & ordering information

Vacuum Range:	0-82Kpa adjustable
Vacuum Display:	dual 'LED' digital display
Vacuum Display Error:	Lower than $\pm 5$ Kpa
Vacuum Stability:	Lower than 5Kpa
Max Free-air Capacity:	8L/min
Time Range:	0.1~99.9 sec
Time Display:	3-digit 'LED' display
Time Set:	0.1 sec
Dimension(mm):	180H x 94W x 260D
Capacity of Collection Bottle:	2000mL
Voltage:	CA220V $\pm 10$ 50~60Hz