

ARCS-100



ARCS-100 Small Shaker

Adjustable speed from 0 to 400 rpm.

“Start” or “Stop” operation can be controlled by front panel function keys or “Start” controlled by connection key while the “Stop” operation controlled by timer flat forms function.

Many kinds of shakers can be exchanged. Can be used inside of the temperature chamber with the max temperature of 80°C.

Type : ARCS-100

Shaking pattern: Reciprocating

Range of varibration : 20mm

Shpg.Wt : 5kg

Speed range : 0-400rpm

Dimensions: 300*360*175