

AS2060B

The AS2060B ultrasonic cleaner characteristics:

Ultrasonic frequency of 40KHz;
Power and smooth;
The strength of two tranches of power conditioning;
Mechanical timer;
Engineering plastic chassis, high-strength aluminum alloy and stainless-steel housing (specify) Optional;
Attached to the grid and noise reduction cover.



Model:	AS2060B
Dimensions (mm):	200W x 185L x 170H
Slot dimensions (mm):	150W x 135L x 100D
Volume (L):	2
Power (W):	60
Power adjustment (%):	50/100
Ultrasonic frequency (KHz):	40
Degassing adjustment:	no
Timing (min):	5 ~ 30 or continuous
Bleeder:	
Housing:	plastic