

AS3120A

AS3120A ultrasonic cleaning device characteristics

Ultrasonic frequency of 40KHz;

The degassing effect of regulatory function, especially for the different viscosity of the solvent degassing;

Electronic timer accurate and durable;

Attached to the grid and noise reduction cover.



Model:	AS3120A
Dimensions (mm):	265W x 160L x 240H
Slot dimensions (mm):	235W x 135L x 100D
Volume (L):	3
Power (W):	120
Power adjustment (%):	30 to 100
Ultrasonic frequency (KHz):	40
Degassing adjustment:	
Timing (min.) :	1-18 or continuous
Bleeder:”>	