

## AS5150B/BD

40KHz/60KHz or other special frequency two steps of conversion is more suitable for different purposes (with the letter D models);

40KHz ultrasonic frequencies can provide strong power;  
60KHz reflect the low noise and fine cleaning;

The degassing effect of the regulatory function is particularly suitable for the different viscosity of the solvent degassing;  
Incidental grid and noise reduction cover



Model:	AS5150B/BD
Dimensions (mm):	320W x 175L x 290H
Slot dimensions (mm):	300W x 150L x 150D
Volume (L):	6
Power (W):	180
Power adjustment (%) :	1-100
Ultrasonic frequency (KHz):	40 (AS5150B)
Ultrasonic frequency (KHz):	40/60 dual-band (AS5150BD)
Timing (min) :	1-999
Bleeder:	
Housing:	aluminum