

AT-330 Column Heater

The unified design makes this column heater's arrangement flexible with the horizontal or vertical setting.

Accredited temperature control and display: Two four-digit LED display can make Temperature control accuracy to $+0.1^{\circ}\text{C}$.

Fast and easy installation of the column: Magnetic elements create an easy opening and closing of the column heater's cover.

Easy view: Crystalline window allows easy watch inside the HTLC column.

Safety free-leaking design: A box placed under heating chamber can prevent liquid from leaking.

Flexible column placement: Provide stainless clamp holder inside heating chamber to make it easy for installing the column under heating.

Smooth pipe: Provide a smooth connection between input and output ports to guarantee that liquid passes through.

Stable column heater base: Provide additional large

size base that can be mounted with column heater once required to stabilize the column heater When placing it vertically.

Convenient sample inserting valve component kits: Sample-inserting valve can be fixed at any height of the column heater.



.Type: AT-330

Temperature: 110°C

Temperature setting range: 0.1°C

Temperature Accuracy: $\pm 0.1^{\circ}\text{C}$

Temperature stability: $\pm 0.1^{\circ}\text{C}$

Caving Size: $360\text{H} * 45\text{W} * 33\text{D}$

Power: 125W

Column installed: 1~2