

BETS-010 Super Thin Orbital Decoloring Shaker



The LCD timing speed display is easy to read and easy to set up. The Brushless Direct Current Motor is efficient and made for long-lasting. It also prevents from electromagnetic interference. Various voltage available and works between 80-240 Voltage. The super engineering ABS plastic molding and platform is elegant and made for long-lasting. It is also environmental-friendly and made without any spray paint. The unique design of the on/off button, the speed knob, the timing knob, and the power cord is fenced with plastic that is resistant from liquids, therefore avoiding any kind of short-circuit power off.

Voltage:	80/240V
50/60HZ	
Power input:	
40W	
Rotation type:	Orbits
horizontally	
Rotation radius:	
15mm	
Frequency:	100 – 350 rpm
Speed and timing display:	LCD
Timing range:	99.59 hours
Platform dimensions:	330 x 330 mm
Overall dimensions:	420 x 370 x 110 mm