





Cell-06							
Picture and Name	Outside Dim H*W*D(mm)	Path Length	Inside Width	Volume	Product Number		
		(mm)	(mm)	(ml)	ES Quartz Glass	IR Quartz Glass	Optical Glass
					190nm-2500nm	250nm-3500nm	340nm-2500nm
 Micro cell with frosted walls and with lid	45*12.5*7.5	5	2	0.35	Q-243	I-243	G-243
	45*12.5*12.5	10	2	0.7	Q-244	I-244	G-244
	45*12.5*22.5	20	2	1.4	Q-245	I-245	G-245
	45*12.5*32.5	30	2	2.1	Q-246	I-246	G-246
	45*12.5*42.5	40	2	2.8	Q-247	I-247	G-247
	45*12.5*52.5	50	2	3.5	Q-248	I-248	G-248
	45*12.5*102.5	100	2	7	Q-249	I-249	G-249
 Micro cell with frosted walls and with lid	45*12.5*7.5	5	3	0.5	Q-253	I-253	G-253
	45*12.5*12.5	10	3	1	Q-254	I-254	G-254
	45*12.5*22.5	20	3	2.1	Q-255	I-255	G-255
	45*12.5*32.5	30	3	3.15	Q-256	I-256	G-256
	45*12.5*42.5	40	3	4.2	Q-257	I-257	G-257
	45*12.5*52.5	50	3	5.25	Q-258	I-258	G-258
	45*12.5*102.5	100	3	10.5	Q-259	I-259	G-259
 Semi-Micro cell with frosted walls and with lid	45*12.5*7.5	5	4	0.7	Q-263	I-263	G-263
	45*12.5*12.5	10	4	1.4	Q-264	I-264	G-264
	45 x12.5*22.5	20	4	2.8	Q-265	I-265	G-265
	45*12.5*32.5	30	4	4.2	Q-266	I-266	G-266
	45x12.5*42.5	40	4	5.6	Q-267	I-267	G-267
	45*12.5*52.5	50	4	7	Q-268	I-268	G-268
	45*12.5*102.5	100	4	14	Q-269	I-269	G-269
 Semi-Micro cell with frosted walls and with lid	45*12.5*7.5	5	5	0.85	Q-273	I-273	G-273
	45*12.5*12.5	10	5	1.7	Q-274	I-274	G-274
	45*12.5*22.5	20	5	3.5	Q-275	I-275	G-275
	45*12.5*32.5	30	5	5.2	Q-276	I-276	G-276
	45*12.5*42.5	40	5	7	Q-277	I-277	G-277
	45*12.5*52.5	50	5	8.75	Q-278	I-278	G-278
	45*12.5*102.5	100	5	17.5	Q-279	I-279	G-279