

DKS-11 Series Water Bath



It adopts stainless steel core and movable top lid. It is precise in temperature control and is well equipped with digital temperature display and chasing alert system.

Specification:

Model	DKS-11	DKS-12	DKS-14	DKS-16	DKS-24	DKS-26	DKS-28
Power supply	220V 50Hz						
Consumption power	300W	500W	1000W	1500W	1000W	1500W	2000W
Temperature control range	RT+10 ~99°C						
Temperature fluctuation	± 1°C						
Follow-up alarm	+2°C						
Dimensions of working chamber(mm)	150*150*110	330*150*110	600*150*110	900*150*110	300*300*110	450*300*110	600*300*810
Overall dimension(m)	300*192*195	450*192*195	750*192*195	1050*192*195	450*345*195	600*345*195	750*345*195