

DN-10 Series Education Digital Microscope



- Digital microscope of integration, 1.3 M pixel digital camera system, clear image and real color reproduction.
- With professional image software for digital microscope, support multilingual operating system, including preview, taking photos or video, image processing and image measurement.
- Friendly user interface, easy operation, the database of colour schemes also can be backuped.
- Widely used in biology, histology, bacteriology etc.

Specification :

			DN-10A	DN-10B		
Viewing Head	Microscope	Monocular head, Inclined at 45°, 360°rotatable		●		
		Binocular head, Inclined at 45°, 360°rotatable			●	
	Camera System	Valid Pixel	1280×1024(1.3M Pixel)		●	●
		Output Mode	USB2.0		●	●
		Operation System	WINDOWS 2000/XP/VISTA 32 Bit or 64 Bit		●	●
		Software	Scopelimage 9.0		●	●
			Scopelimage Advanced		●	●
		Rang of Viewing Field	> 75%		●	●
Objective	Achromatic Objective (DIN) 4×、 10×、 40×、 100×			●	●	
Eyepiece	WF10×/ 18			●	●	
Stage	Double Layers Mechanical Stage 115×125mm/70×20mm			●	●	
Condenser	Abbe NA1.25 Condenser with iris diaphragm (Φ2-Φ28mm)and frosted filter			●	●	
Focusing	Coaxial coarse and fine adjustment			●	●	
Illumination	1W Powerful LED Illumination, brightness adjustable			●	●	

Note:● Standard Outfit, ○ Optional