

DN-107T Digital Microscope



- Integrated digital microscope, resolution 1280×1024 (1.3 pixel)
- USB2.0 output, plug and unplug freely
- Support WINDOWS 2000/XP/VISTA system, Equipped with powerful applied software
- Economic features, widely used in academic demonstration and health system

		DN-107T		
Viewing Head	Sliding Binocular Head, Inclined at 45°, Interpupillary 55-75mm		●	
	Camera System	Valid Pixel	1280×1024 (1.3M Pixel)	●
		Output Mode	USB2.0	●
		Operation System	WINDOWS 2000/XP/VISTA	●
		Software	ScopeImage 9.0	●
		Range of Viewing Field	>75%	●
Eyepiece	WF10×/ 18		●	
	P16×/11		○	
Nosepiece	Toward Quadruple Nosepiece		●	
Objective	Achromatic Objective 4×、 10×、 40×、 100×		●	
Stage	Double Layers Mechanical Stage 140×140mm/75×50mm		●	
Condenser	Abbe NA1.20 with Iris Diaphragm & Fliter		●	
Focusing	Coaxial Coarse & Fine Adjustment System, Range 28mm, Fine Division 0.002mm		●	
Illumination	Halogen Lamp 6V/20W, Brightness Adjustable		●	
Phase Contrast Kit			○	
Dark Field Attachment			○	
Polarization Attachment			○	

Specification :

Note:● Standard Outfit, ○ Optional