

Db-Va Digital Display Thermostatic Hotplate



Basic Info

- Model NO.: DB-VA
- Digital Display: Stainless Steel Plate
- Specification: 1 unit in 1 box
- HS Code: 8419899090

- Type: Heating Type
- Trademark: OEM/Unbranded
- Origin: China Mainland

Product Description

DB-VA Thermostatic Stainless-Steel Laboratory Hotplate

Laboratory Hotplate

1. Stainless Steel heating plate, hi-power, uniform heating, LED
2. CE, ISO13485/14001/9001
3. Factory

The factory mainly produces various kinds of teaching and studying instruments medical and epidemic preventive instruments, serial analyzing and experimental instruments and so on over one hundred of products.

Specifications of DB-VA Thermostatic Stainless-Steel Laboratory Hotplate	
Power Supply	AC 220V/110V 50Hz/60Hz
Heating power	2000W
Temperature range	Room temperature to 300 Celsius
Overall size	L40 W68 H12cm
Working area	L40W60
Digital display	Yes