

GZ-GZ-A



GZ(GZ-A) High Power Laboratory Overhead Stirrer

High Power Laboratory Overhead Stirrer, Hi-power, 20L capacity, 2-step regulate, durable

Spec:

GZ(GZ-A) High Power Electric Stirrer/Overhead Stirrer	
Power Supply	AC 220V 50Hz/60Hz
Speed	60-500rpm(stepI) ,240-2000rpm(step II)
Motor power	200W
Max. stirring capacity	20Liter
Maximum moment of stirring shaft	1850gcm
Applicable medium viscosity	0-10000mpas
Large capacity, easy to use, cost efficient, widely applicable.	