


HD-1 THICKNESS TESTER



HD-1 THICKNESS TESTER

HD series of instruments are used for detecting the thickness of tablet and capsule.

Applicable standard

Corporate standard(Thickness tester) Q/12XQ0194-2010

Features:

Used for measuring the thickness of peak and wall of the capsule.

High-accuracy test with the most precise, high-quality homemade centigrade meter.

Specifications:

Test range: (0~5)mm

Test accuracy: ± 0.01 mm

Type of capsule: 00,0,1,2,3,4 type capsule

Dimension: (200×80×50)mm³