

HN-AT5A Low Speed Large Capacity



Bench top, large capacity, horizontal rotor, brush less motor, multi-capacity, sound structure, durable a large range of applicability's.

HN-AT5A Low Speed Large Capacity Centrifuge

Technical data:

Speed: start-5000 rpm (touch key control)
Capacity: 250ml×4, 50ml×8, 15ml×24 (3 rotors)
Max RCF: 4300 xg
Timer: 0-99min
Power Supply: 220V/110V 50Hz/60hz
Power: 320W