

# HZQ-3112-SERIES



## Microprocessor controller(with timing function)

### Features:

- \* Large LCD screen to display more data at same time, menu-driven interface easy to understand.
- \* 304 stainless steel platform, easy to remove and clean.
- \* Automatically saves the date in case of power failure, and run the previous programs after power on Microprocessor controller with timing function and slow starting function to prevent the liquids spilling, also avoid equipment unexpected start and stop Unique design of speed monitoring circuit, it will stop shaking when the speed is too higher too low.
- \* DC motor with low heat.
- \* Three eccentric shaft balancing drive ensures that the sample on platform shaking with same speed.
- \* Programmable controller is option.

Spec:

Model	HZQ-3112 (Double-desk)	HZQ-3111 (Monolayer)	HZQ-3222 (Double-desk)	HZQ-3221 (Monolayer)
Display Resolution	LCD			
Rotating Speed	20-300rpm			
Speed Accuracy	+/-1r/min			
Amplitude	35mm/50mm			
Timing Range	0~500h			
Standard Configuration	500ml*60PCS		500m*60PCS	
Maximum Configuration	250ml*48pcs or 500m*30pcs or 1000m*20pcs or 2000ml(12pcs (every desk)		250ml*63pcs or 500ml*35pcs or 1000ml*35pcs or 2000ml*16pcs (every desk)	
Platform Size	800mm*600mm		920mm*640mm	
Platform Shelves	2	1	2	1
Exterior Dimension	900mm*680mm*680/mm		1060mm*850mm*680/mm	
Power Consumption	400W			
Electrical Requirement	AC220V			