

DC Digital Display Constant Speed Power-Driven Mixer

Series DC digital display constant speed electric mixer has reasonable design and novelty structure. It use DC brushless motor and step less speed, and is equipped with voltage stabilizing system, so the display of rotation speed is stable. It has such features as large power, little noise, safe operation, and soon, which had provided one indispensable basic equipment for mixing various liquid-phase chemical reactions at labs of colleges and schools, scientific research unit and chemical plant.



JJ-4 Six-link
Synchronization
Electric Mixer



JJ-6 Digital Display
DC Electric
Mixer



JJ-2 Force
Increasing



JJ-1 Precision
Electric Mixer



GZ (GZ-A) High Power
Laboratory Overhead
Stirrer

| Type | Description | Power of Motor | Timing Range | Rotation Speed | Speed Digital Display | Weight |
|----------|---|--------------------------------|--------------|----------------|-----------------------|--------|
| JJ-1JJ-2 | Force Increasing /Electric Mixer | 25w,40w,60w,80w,100w,150w,200w | 0-120min | 0-3000rpm | No | 10kg |
| JJ-3JJ-4 | Six-link Electric Mixer | 6*100w | 0-120min | 0-1800rpm | Yes | 32kg |
| JJ-6 | Digital Display DC Constant Speed Mixer | 100W,300W | 0-120min | 0-2500rpm | Yes | 12kg |

| GZ (GZ-A) High Power Electric Stirrer/Over head Stirrer | |
|---|--|
| Power Supply | AC220V50Hz/60Hz |
| Speed | 60-500rpm(StepI),240-2000rpm(StepII) |
| Motor power | 200W |
| Max. stirring capacity | 20Liter |
| Maximum moment of stirring shaft: | 1850gcm |
| Applicable medium viscosity | 0-10000mpas |
| Display | GZ-A is digital, GZ is without display |
| Large capacity, easy to use, cost efficient, widely applicable. | |