

# JY02S UV Transilluminator



## Technical specifications

Transmission size ( W×L ) : 200×150mm

Source ultraviolet wavelength: 302 ( nm )

Reflectance UV wavelength: 254 、 365 ( nm )

Transmission UV lamp power: 8 ( W )

Reflectance UV lamp power: 254nm ( 11W ) 、 365nm ( 11W )

Volume ( L×W×H ) : 435×295×490 ( mm )

Weight: 16.9kg

## Characteristic feature

- No need for darkroom, can be use for all-weather.
- With the operation of export for cutting the gel.
- Using the UV filter glass, quartz ultraviolet lamp, long service life, light uniform.