

KDM-SERIES



KDM Thermostatic Sheath

Model	KDM	KDM	KDM	KDM	KDM	KDM
Capacity	250ml	500ml	1000ml	2000ml	3000ml	5000ml
Heat Power	200W	400W	500W	800W	1000W	1500W
Power	AC 220V					

KDM-A Digital Display Thermostatic Sheath

Model	KDM-A	KDM-A	KDM-A	KDM-A	KDM-A
Capacity	250ml	500ml	1000ml	2000ml	3000ml
Heat Power	200W	400W	600W	800W	1000W
Temp. Range	RT-300°C				
Power	AC 220V				