

KJ201C Microplate Shaker



Features:

- © Elastic silicon plate, easy to clean and wear-resisting. It can fit various kinds of cell plates. 6pcs of cell culture plates can be mixed at one time. High performance & low noise AC motor to make rapid and thorough mixing.
- © Adjustable knob design can ensure step-less speed adjustment and constant run within 60 min.

Technical Data:

Operation method: Circling
Circling diameter: 9mm
Voltage: AC 220V \pm 10% (50/60Hz)
Power rating: 10 W
Oscillating Frequency: 200-2800rpm
Timing range: 0-15min
Operate mode: Timing / Continuous
Max. loading: 6 pcs microplate
Working plate size (6 plates): 278 \times 270mm
(1 plate): 128 \times 86mm
Dimensions: 310 \times 280 \times 225mm
Net. Weight: 6.2Kg