

KJMR-V 3D Rotation Mixer

The kind of Mixer with three-dimensional motion provides gentle and thorough mixing. It is ideal for sample mixing in centrifuge tubes, as well as resuspension of large, fragile DNA FRAGMENTS. The KJMR-V accommodates most size tube, centrifuge tube, microwell plate, petri dishes and tissue culture plate and cell culture bottle. This equipment is suitable to work in the incubator or cold room.

Main Technical Data:

Voltage: 220V±10%
(50/60Hz)

Power Rating: 6W

Type of Movement: 3D

Angle: 20°

Maximum Load Capacity:

0.8Kg External Dimension:
240*170*150mm

