

# LCD-600



LCD-600

## FEATURES

- 1.The image can be easily viewed by many people on monitor at the same time.
- 2.Improve the working condition also can soften the tiredness caused by long time observation.
- 3.Integrated structure,easily operation,small size and can save space.
- 4.The magnification can be zoomed non-step.
- 5.There is “+” for purpose of positioning and measurement.
- 6.Easily equipment,it can be work properly by put it on a stable place.

## OPTICAL DATA

Model	Objective Magnification	Working Distance	Optical Magnification ( 1/3CCD with 8"inch LCD screen)
LCD-600	0.65×-6.0×	100 mm	11×-102×

## SPECIFICATION LIST

Components	Specification	Model
------------	---------------	-------

		LCD-600-Y	LCD-550-W	LCD-450
LCD screen	8"英寸 8"inch	●	●	●
	10.4"英寸10.4"inch	○	○	○
Reduce lens	0.5×			●
	0.45×	●	●	
Objective	0.75×-6.0×	●		
	0.65×-5.5×		●	
	0.7×-4.5×			●
3-D viewing lens	Manual, 360° viewing	○	○	○
Laser locator	Red	○	○	○
Resolution	1/3" 420 线 1/3" 420 line	○	○	●
	1/3" 520 线 1/3" 520 line	○	○	○
	1/3" 580 线 1/3" 580 line	●	●	○
Incident light	Halogen lamp 12V/10W		●	●
	96T LED ring light	○	○	○
Transmitted light	Fluorescent Lamp 5W		○	○
	Halogen Lamp 12V/10w		●	●