

LX-100MiniCentrifuge

Characteristic:

- * It is specially designed by Molecular Biologist, it can be held in palm to operate.
- * It is used for the operation with normal low-speed, and it can reduce the usage frequency of high-speed centrifuge;
- * It can be used as the special centrifuge for isotope and it can reduce the contaminated rate of high-speed centrifuge;
- * The transparent arch cover can be easy to observe the operation;
- * The press switch is very convenient to use, it will begin to run after close the cover, and when you open the cover, it will stop automatically;
- * It is specially applicable for the low-speed centrifugal operation which collect the liquor around EPPENDORF TUBE to the bottom;
- * Built-in temperature relay, it will turn off the power supply after long time operation to protect the usage life of motor;

Specification:

Rotor Capacity	1.5ml*6	Rotor:	1.5ml 1pc
	0.5ml*6		0.5ml tube 10pcs
	0.2ml*6		0.2ml tube 10pcs
Dimensions	300*335*275mm(l*w*h)	Weight	0.5KG
Power	220V 50Hz 30W (110V 60Hz Option)	Max. Speed	3000rpm

