

Quadruple Beam Balance

This product is widely used for weighing in industrial enterprises (such as chemical, medical and machine-building enterprises), agricultural applications and schools. Outer packing:
Carton dimensions: 720 x 450 x 380mm
Quantity: 4PCS



Model	MB311
Capacity	311g
Readability	0.01g
Beam Calibration	FrontBeam1*0.01g 2" Beam10*1g 3" Beam100*10g 4"Beam200*100g
Damping	Magnetic
PanSize	100mm(4")diameter Stainless steel
Density platform	75mm (3")diameter Adjustable height
Overall size(w*h*d)	370*300*137mm(14.5**11.8**5.4")
Packing dimensions(L*W*H)	720*450*380(4pcs)
Gross weight	10KGS(4pcs)