

## MH-2810D



2810D-type multi-use characteristics of incubators

Increase in the 2800A on the basis of the timing functions and timer time can be set between 1-999 minutes  
Regular time set by the buttons on the front panel by LED digital display.  
Can accommodate larger heating block.

Part Number: 2810D

Temperature Display: dual four-digit display

Constant temperature range (°C): room temperature +5-200

The thermostatic stable degrees (° C):  $\pm 0.1$

Thermostat accuracy (° C):  $\pm 0.15$

Upper limit temperature protection (°C): 200

Temperature deviations protection (°C): set temperature +5

Timing (min) : 1-999

The incubator Size (mm): 215W x 51H x 135D

To accommodate the heating block number: 2 (standard) / (double)

Heating power (W / V / Hz): 750/220/50

Dimensions (mm): 320W x 92H x 270D