

## MH-5800



Constant temperature heating coil using a high-precision PID temperature controllers, temperature control is fast and smooth, and advanced user self-tuning function. The heating coil of the thermostat can provide the platform of a large area of constant temperature. It can be widely used in biochemical, environmental, pharmaceutical and other fields.

### Features

- \* PID temperature controllers, high precision temperature control, and temperature stability.
- \* Pairs of four LED digital display to set the temperature value (SV) and the actual temperature value (PV).
- \* Has over-temperature alarm power-off function and safe operation.
- \* “Pu-b \* with user operable temperature deviation correction.
- \* Heating coil disk is made of stainless steel.
- \* Ultra-thin design, beautiful, user-friendly.

Model: MH-5800 multi-purpose thermostatic plate

Temperature range Temperature range: room temperature for more than 3 ° C ~ 80 ° C

Heating power: 240W

Temperature accuracy: better than  $\pm 0.5$  °C

Dimensions (LxWxH) (mm): 515 x 325 x 45

Workspace size (LxW) (mm): 300 x 400

Power supply (V / Hz): 220/50

Mains fuse tube (A): 3