

MiniT-1/MiniT-3 Biological Indicator Incubator

The MiniT-1/MiniT-3 Biological Indicator Incubator is designed for temperature incubation of the EO Biological Indicators. It provides easy verification of incubation temperature.

Features:

1. LCD display, with 12V power for ease of portability
2. The incubator has many wells for placement of the biological indicator vials and the wells can be custom-made. User can set the incubating temperature.
3. Materials: Aluminium block
4. Well Diameter: According to custom-making

Specification:

Model	MiniT-1	MiniT-3
Temperature Control Range	RT+5°C~80°C	RT+5°C~100°C
Timing Range	1min~99h59min	
Temperature Control Accuracy	±0.5°C	
Display Accuracy	0.1°C	
Heating Time (20°C to 57°C)	≤6min	
Temp. Control Accuracy @ 57°C	0.35°C	
Dimensions (LxWxH) mm	110x150x80	120x152x112
Net Weight	0.5kg	0.85kg

Block (Option)

A-MiniT-1	15xΦ10mm
B-MiniT-1	12xΦ12mm
C-MiniT-1	6xΦ15mm
D-MiniT-1	6xΦ16mm
E-MiniT-1	18xΦ9mm

