

MTN-2800D



TYPE: MTN-2800D-12/24/40
TEMPERATURE RANGE AMBIENT: +5°C~180°C
TEMPERATURE UNIFORMITY: $\pm 0.1^{\circ}\text{C}$
POWER-OFF TEMP. : 185°C
DEVIATION TEMP. : SET TEMP+5°C
TIMER RANGE: 10MIN~99HR59MIN
HEATED ZONE(MM): 95WX70HX75L(STANDARD SIZE)
150WX70HX95L(double size)
NUMBER OF HOLES: 6/12/20(STANDARD SIZE)
12/24/40(DOUBLE SIZE)
HEATING POWER: 350W
MAX BLOCK : 1(STANDARD SIZE)
Dimensions(mm) : 240WX130HX215D
Power(VAC,Hz): 220V/50Hz