


MTN-2810D



TYPE: MTN-2810D-8/12

TEMPERATURE RANGE AMBIENT : +5°C~300°C

TEMPERATURE UNIFORMITY : $\pm 0.1^{\circ}\text{C}$

POWER-OFF TEMP. : Choose as your requirement

DEVIATION TEMP. : Choose as your requirement

TIMER RANGE : 10MIN~99HR59MIN

NUMBLE OF HOLES : 8/12

HEATING POWER : 1350W

HEATED ZONE(MM) : 210WX95HX150L(STANDARD SIZE)

Dimensions(mm) : 268WX230LX165H

Temperature : 5°C~35°C

Power(VAC,Hz) : 220V/50Hz