


# N-117M Biological Microscope



New machining facility and advanced alignment technology

- Complete functions and outfit of binocular microscope
- Comfortable operation with updated and ergonomic design

## Features:

	Comfortable observation with adjustable eye point height and interpupillary distance
--	--

## Specification :

		N-117M
Viewing Head	Compensation Free Binocular Head, Inclined at 30°, Interpupillary Distance 48-75mm	●
	Compensation Free Trinocular Head, Inclined at 30°, Interpupillary Distance 55-75mm	○
Eyepiece	WF10×/18	●
	P16×/11	○
Objective	Achromatic Objective 4×、 10×、 40×、 100×	●
Nose piece	Backward Quadruple Nose piece	●
Stage	Double Layer Mechanical Stage 132×142mm/ 75×40mm	●
Focusing	Coaxial Coarse & Fine Adjustment, Moving Range 24mm, Fine Division 0.004mm	●
Condenser	Abbe NA 1.20 with Iris Diaphragm & Filter	●
	Abbe NA 1.25 with Iris Diaphragm & Filter	○
Illumination	S-LED, Brightness Adjustable	●
	Halogen Lamp 6V/ 20W, Brightness Adjustable	○
Phase Contrast Kit		○
Dark Field Attachment		○
Polarization Attachment		○

Note: ● Standard Outfit, ○ Optional