

N-200M Biological Microscope



With excellent optical system and reliable structure, it makes the image much sharper, and the operation easier. With wide range of accessories [e.g. dark field, phase contrast, epi-fluorescent] and advances to expand the application scope of the N-200M, the versatility of the microscope provides More Microscopy for you.

Features:

| | |
|---|---|
|  | Big Operating Space with Backward Nosepiece, Sharp Image with High Quality Objectives |
|  | Easier Comparing Test with Dual Slides Finger of Mechanical Stage |
|  | Comfortable Handling with Precise Focusing System |

Specification:

| | | N-200M |
|----------------------------|--|--------|
| Viewing Head | Compensation Free Binocular Head, Inclined at 30° | ● |
| | Compensation Free Trinocular Head, Inclined at 30° | ○ |
| | Monocular Head, Inclined at 30° | ○ |
| Eyepiece | Wide Field Eyepiece WF10×/ 18 | ● |
| | Wide Field Eyepiece WF16×/ 11 | ○ |
| Objective | Semi-plan Achromatic Objective 4×, 10×, 40×, 100× | ● |
| | Plan Achromatic Objective 4×, 10×, 20×, 40×, 60×, 100× | ○ |
| Nosepiece | Backward Quadruple Nosepiece | ● |
| | Quintuple Nosepiece | ○ |
| Stage | Double Layers Mechanical Stage 140×140mm/ 75×50mm | ● |
| Condenser | NA 1.25 Abbe Condenser with Iris Diaphragm and Filter | ● |
| Illumination | S-LED Illumination, Brightness Adjustable | ● |
| | Halogen Lamp 6V/ 20W, Brightness Adjustable | ○ |
| | Kohler Illumination | ○ |
| Filter | Blue Filter | ○ |
| | Green Filter | ● |
| Polarization Set | Simple Polarization Set | ○ |
| Phase Contrast Kit | Simple Phase Contrast Kit | ○ |
| | Sliding Phase Contrast Kit | ○ |
| | Turret Phase Contrast Kit | ○ |
| Dark Field Attachment | Dark Field Condenser (Dry) NA 0.9 | ○ |
| | Dark Field Condenser (Oil) NA 1.25-1.36 | ○ |
| Epi-fluorescent Attachment | | ○ |
| Pointer | Green LED Pointer System | ○ |

Note: ● Standard Outfit, ○ Optional