

QB-9001



QB-9001 Microporous Quick Shaker

Breaking through convention with first-class originality and updated structure. Streamline main frames formed by enhanced type engineering plastics, up to 4 to 6 mini-plate quick shaking, non-skid natural rubber shaking plate, and timing and continuity all together bring you an optimum platform for experiments.

Spec:

Electricity: 220V

Capacity: 30W

Speed: Step-less speed regulation

Frequency: 0 – 1350 rpm

Timing: 0~120 min / continuous

Loading: Four 96-pore mini plate

Overall dimensions: 350×250×150mm