

R1002B-10L Rotary Evaporator



Overview

10L Rotary Evaporator is the pilot-scale unit to offer vacuum tight, uniform speed of rotation, and the constant temperature of the heating bath, which is robust, durable and user-friendly.

- – Capacity: 10L;
- – Teflon seal, ultimate vacuum <2 mbar (mechanical sealing is option);
- – Heating bath, PID temperature control range: RT ~ 99°C (180 °C, Oil Bath);
- – Motorized lift mechanism;
- – Quality borosilicate glass 3.3 (BG3.3);
- – Vertical dual coil condenser;
- – Rotary speed: 10 – 110 rpm, LED display;
- – Vacuum meter equipped;
- – Teflon discharge valve for clean discharge.