

SB-3200DT Ultrasonic Cleaner



Features:

- Digital display and setting ultrasonic cleaning time
- Display working countdown time
- Display actual temperature of the cleaning tank
- Temperature setting range: Ambient temperature-80°C
- Timing range: 1-99 minutes adjustable
- Equipped with a special stainless-steel basket, reduce sound cover
- The shell and cover are made of stainless steel
- The equipment controlled by Single-chip controller
- The latest working parameters can be stored in memory after power-off

Main Technical parameters:

Model	Inner Tank	Volume	Frequency	Power	Heating power	Adjustable Temperature	Adjustable Time	Drainage	Rack
	L/W/H (mm)	L	KHZ	W	W	°C	min		
SB-100DT	300·150·100	4.5	40	100	150	Ambient-80	1-99	/	Y
SB-120DT	240·140·100	3	40	120	150	Ambient-80	1-99	/	Y
SB-120DT	240·140·150	5	40	120	150	Ambient-80	1-99	/	Y
SB-3200DT	300·155·150	6	40	180	180	Ambient-80	1-99	/	Y
SB-5200DT	300·240·150	10	40	200	320	Ambient-80	1-99	Y	Y
SB-5200DT	300·240·150	10	40	250	320	Ambient-80	1-99	Y	Y
SB-5200DT	300·240·150	10	40	300	320	Ambient-80	1-99	Y	Y
SB-4200DT	300·240·200	14.4	40	400	320	Ambient-80	1-99	Y	Y
SB25-12DT	500·300·150	22.5	40	500	800	Ambient-80	1-99	Y	Y
SB25-12DT	500·300·150	22.5	40	600	800	Ambient-80	1-99	Y	Y
SB-800DT	500·300·200	30	40	800	800	Ambient-80	1-99	Y	Y
SB-1000DT	600·300·300	54	28	1000	2800	Ambient-80	1-99	Y	/
SB-1200DT	600·400·300	70	28	1200	3200	Ambient-80	1-99	Y	Y
SB-1500DT	600·500·300	90	28	1500	4000	Ambient-80	1-99	Y	/
SB-1800DT	600·600·300	108	28	1800	4000	Ambient-80	1-99	Y	/
SB-2000DT	650·600·300	117	28	2000	5000	Ambient-80	1-99	Y	/
SB-3600DT	750·600·300	135	28	3600	6000	Ambient-80	1-99	Y	/
SB-6000DT	1200·650·350	273	28	6000	8000	Ambient-80	1-99	Y	/
SB-8000DT	1800·650·350	410	28	8000	10000	Ambient-80	1-99	Y	/