

SHZ-82 & SHZ-82A



Shaking Water Bath

Durable, 2000W, SUS inner, timer, 12kg of the load.

The thermostatic shaking water bath is also called as thermostatic water bath swing bed. The temperature of this serial product adopts integrated operational amplifier and LED display, which has high accuracy to control temperature, easy adjustment, correct and direct values to be displayed, superior and reliable functions. The incubator of the water bath constant temperature oscillator adopts the imported stainless steel and advanced process for production have the characteristics of stable oscillation, reliable quality, long time, etc. Gripper adopts advanced stainless steel spring Multi-purpose fixture, convenient to use.

This serial product is suitable for production and scientific research departments such as environmental protection, medical treatment, teaching and studying, sanitary and epidemic prevention, pharmaceutical inspection, zoology and botany, ocean science, foodstuff engineering etc. It is a special constant temperature equipment with oscillating for BOD measurement of water body analysis; For cultivation and preservation of bacteria, virus, mold, and microbe; For plantation of plants and breeding test.

Specifications of SHZ-82/SHZ-82A Laboratory Thermostatic Shaking Water Bath

Power Supply:	AC 220V/110V 50Hz/60Hz
Temperature control range:	Room temperature to 100 celsius
Temperature fluctuation:	1°C
Amplitude:	20mm (Convoluting oscillation)
Heating power:	2000W
Timer:	0-120min
Working load:	12kg
Oscillation tray:	410x310mm
Capacity:	49x39x17cm
Overall size:	L72xW48xH52cm
Digital display:	LED digital display temp. control
Oscillation frequency:	SHZ-82:0-280times/min stepless regulation (Orbital)
	SHZ-82A: 0-280times/min stepless regulation (Rotary)
Easy to use, quality stainless steel material, good reliability.	