


## TDL-24/TDL-24C



### Product Description

It is possessed of these special feature following: It is beautiful in appearance, large capacity but small in bulk and has a complete function and a stable performance, you can choose the speed that you need at the set range, and it can keep balance automatically, in addition, it has a little rate of increase in temperature and can be used efficiently and widely, It is a dependable instrument for qualitative analysis of serum, plasma, and immunity in hospitals or chemical and biochemical laboratories.

Easy to operate, with large capacity, steady performance. Applicable for blood serum and plasma separation in laboratories.

TDL-24 is analog display and TDL-24C is digital display.

### Specification:

Technical data:

Power Supply:	AC 220V 50Hz
Speed:	0-4000 rpm
Timer range:	0-60min
Max. rxf:	2325xg
Capacity:	10mlX24holes or 15mlX24holes