

TDL-60B/TDL-50B/TXL-4.7 LowSpeed Centrifuge

The machine's frame is made of high strength plastic, so it's weight is light, the noise is low when running, it is controlled by microcomputer. The motor has a wider range of speed & runs precisely. It can be operated simply and conveniently. It is safe and reliable because of being equipped with an electronic door lock.

It can be used for qualitative analysis of serum, plasma, evictions in the fields of biology, chemistry, medicine and so on.

Specification:

Timer Range	0~60Min
Noise	<70dB
Power	200V/50Hz(110V optional)

Rotor:

Model	Max Speed	Max. RCF	Rotor Type	Capacity	
No.1 Rotor	6000RPM	5200xg	Angle Rotor	20ml*12	For TDL-60B
No.2 Rotor	5000RPM	3800xg	Swing Bucket Rotor	20ml*8	For TDL-50B
No.3 Rotor	4700RPM	2300xg	Angle Rotor	7ml*12	For TXL-4.7

