

Vortex-BE1 Vortex Mixer



The universal rack is able to fit a variety of centrifuge tubes ranging from 0.5 ml to 50 ml. The powerful and stable mixer is made to meet various experimental requirements with mixing speeds in the 200 – 2800 rev / min.

It is convenient with a three – point switch that can be set to continuous rotary or touch. The double – damping rubber is specially designed to prevent sliding test rig equipment.

Amplitude:	4 mm
Speed range:	200 – 2800 rev / min
Working mode:	Continuous touch and speed regulation
Platform size:	70 mm
Power:	30 W
Dimensions:	145 x 175 x 140 mm
Weight:	2 kg