

## WB-2010A



Micro programmable PID controller reduces the temperature difference inside shaking bath and smoothly raises the temperature.

Two 4-digit “LED” display in the front panel shows the “Set” and “Actual” temperatures.

It can automatically shut the power off when being overheated, making the system safe and reliable.

Temperature offset calibration function provide temperature accuracy control.

Number of the shaking per minute can set.

Operate continuously or set the timer up to 90 minutes.

Special design to prevent additional water drop into flasks.

Provides different kinds of sample holders.

Type: WB-2010A

Heatig Power: 1300W

Capacity: 20L

Tank Dimension (mm): 490W \* 295H \* 150D

Temp. Range: ambient~95°C

Temp. Accuracy:  $\pm 0.2^{\circ}\text{C}$

Outside Dimension: 500W \* 300H \* 300D