

## WRR Visual Melting-Point Apparatus



### Characteristic:

\* To be used to measure the melting point of medicines, dye-stuffs and perfumes and other organic crystal substances of determine their purity. Measurement method is in conformity with the standard of pharmacopeia.

\* Measuring three samples at the same time and calculate average value of initial and final melting points automatically.

### Specification:

Melting-	40~280°C
Minimum digital display:	0.1°C
Linear heating rates Ω	0.5,1.0,1.5,3.0°C/min
Measuring accuracy Ω	<200°C;±0.5°C.>=200°C;±0.8°
Repeatability (standard Deviation ó)	0.3£C
Weight(Net):	12.5Kg
Overall Dimensions:	270mm*310mm*400mm