

## WRS-1B Digital Melting-Point Apparatus

- \* Measurement of melting point is one of the important methods of measurement of purity.
- \* To be used in production of medicines, dye-stuffs, perfumes and other organic crystal substances.
- \* Photo electric test
- \* LCD display.
- \* RS-232 Interface and real time melting curve.

### Specification:

Melting-point measuring Range:	Roomtemp to 300°C
Minimum digital display:	0.1°C
Linear heating rates :	0.2, 0.5, 1.0, 1.5, 2.0, 3.0, 4.0, 5.0°C/min
Measuring accuracy :	<200°C: ±0.5°C , 200°C~300°C: ±0.8°C
Overall Dimensions:	510mm*330mm*200mm

