

## Measure

### WYA-2SABBEDigitalRefractometer

#### Characteristic:

- \* Measurement of refractive index nD of transparent or translucent liquid and solid substances.
- \* Measurement of the Brix of sugar solution.
- \* Visual aim and LCD display
- \* Automatic correction of the effect of temperature on the Brix.
- \* Prism made of hard glass, un easy to be worn.
- \* Printer interface for direct data print-out.

#### Specifications:

Measurement range:	Refractive index nD:1.3000-1.7000
	Dissolved SolidsBrix:0-95%
Accuracy:	Refractive index nD:±0.0002
	Dissolved Solids Brix:±0.1%
Temperature Range(Min.0.1°C):	0-50°C
Weight (net):	10kg
Overall dimensions:	330mm*180mm*380mm



WYA-2SABBEDigital Refractometer