

XSZ-107T Series



Features:

· With precision machining equipment and advanced alignment technics, making improvement on the basis of traditional microscope XSZ-107BN. It is an ideal instrument in biological, histological, pathological field and can be widely used in academic demonstration and research.

Specification:

		XSZ-107T	XSZ-N107T
Viewing Head	Sliding Binocular Head	●	●
	Sliding Trinocular Head	○	○
	Seidentopf Binocular Head	○	○
	Anti-mould	○	●
Eyepiece	WF10X	●	●
	P16X	○	○
Nosepiece	Quadruple Nosepiece	●	●
	Quintuple Nosepiece	○	○
Objective	Achromatic Objective 4X,10X,40X,100X	●	●
Stage	Double Layers Mechanical Stage 140X140mm/ 75X50mm	●	●
Condenser	Abbe NA1.20 with Iris Diaphragm & Filter	●	●
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation, Moving Range 28mm	●	●
Armrest	Black Armrest	○	●
Rubber Cover	Rubber Cover for Focus Adjustment Knob	○	●
Illumination	S-LED Illumination, Brightness Adjustable	●	●
	Halogen Lamp 100V-240V/ 6V20W	○	○
	Plan-concave Mirror	●	●
Package	Styrofoam	●	●
	Carton	●	●
	Wooden Case	○	○

Note: ● Standard Outfit, ○ Optional