

ZTX-3S



ZTX-3S-C2

FEATURES

XTX-6S, ZTX-3S series microscope can provide high resolution and sharp image. With built-in CMOS software, it is very easy for demonstration and analysis. The instruments are widely used in educational field and electronic industry.

OUTFIT

Specification list	Specification	Model	
		XTX-6S-W	ZTX-3S-C2
Viewing Head	Inclined Binocular	●	●
	Head 360° Rotatable	●	●
	Built-in CMOS(1.3 M pixel)	●	●
Eyepiece	WF5×/φ16	○	○
	WF10×/φ20	●	
	WF15×/φ15	○	○
	WF20×/φ10	○	○
	SWF10×/φ123	○	○
	HWF10×/φ22		●
Objective	Conversion type 2×/4×	●	
	Conversion type 1×/2×、 1×/3×	○	
	Zoom magnification 1×-4×		●
Focusing	Coarse Focus Adjustment System	●	●

	Coaxial Coarse & Fine Focus Adjustment System		○
Incident light	Halogen lamp 12V/10W	●	●
	LED Lamp (3.2V20MA)×3	○	
	LED Lamp 3.2V1W	○	○
Transmitted light	Halogen lamp 12V/10W	●	●
	LED Lamp (3.2V60mA)	○	○
Brightness adjustable system		●	●

● **Standard** ○ **Optional**